




Tools for Vaccine Development and Covid-19 Defense



MG Scientific thanks all of you working diligently in support of the fight against COVID-19. Consistent with our commitment to keep you informed about the latest product developments from our manufacturing partners, MG Scientific is pleased to present the following featured products from our friends at Bel-Art, Eppendorf, and Globe Scientific. Click here to see all current [Promotions](#), as well as updated [Covid-19 Information and Solutions](#).



Pipettors



New MOVE-IT Pipettors

Double your pipetting performance. Eppendorf redefined adjustable tip spacing pipettes from scratch according to their premium standards. Move It is the result of customer requests for an efficient and safe solution for synchronous pipetting of a series of samples between different vessel formats, such as between tubes and plates. What is unique about the new Move It? It can do this without any tubing connections between cone and piston-cylinder system. This enhances pipetting performance: Fewer movable, often fragile parts.



The Diamond Apex series of pipettors are simply one of the best lines available today. The magnetic assisted piston ensures consistency of plunger motion to the first and second stops for aspiration, dispense and touch off. The new spring and seal design allows for low plunger force, reducing repetitive fatigue. A corrosion-resistant alloy used in the tip ejector is designed to enhance durability.



Search on-line to see what new and innovative products *MG Scientific* has for your lab or call *MG Customer Service* to see the many different solutions to fit your needs- 800.343.8338