



# Multi-brand Pipettes & Pipet Accessories



The need for accurate, precise, and reproducible liquid handling tools is a common requirement for all laboratories regardless of analytical focus. MG Scientific recognizes, however, that each lab's requirements will be unique to their specific applications. Therefore, MG Scientific offers its customers options. Featured below is a sample of the wide variety of liquid handling tools and accessories that we offer. We encourage you to contact your MG Scientific representative to discuss your specific application needs.

## BrandTech - HandyStep



**HandyStep® Touch & Touch S**  
The first electronic repeating pipettes with touchscreen operation and automatic tip ejection. Available in the basic HandyStep touch model, or HandyStep touch S model with added features.

[Learn More](#)

**CellTreat Reagent Reservoirs**

- 8 and 12-channel
- Holds up to 5mL each
- Compatible with Single and Multi-channel pipettors

[Learn More](#)



## CellTreat - VistaLab



**RADICALLY different**

Fully mechanical volume adjusting Single and Multi-channel Pipettors. audible "click" tip attachment to ensure consistent loading, with easy spring-loaded "flick" tip ejection for low thumb force



[Learn More](#)

## Thermo Scientific - Finnpiptettes, Finntips, & Stand

Use the Thermo Scientific Finnpiptette F1 & F2 Single and Multi-Channel Pipettes for highly demanding applications where robustness and durability are essential. Provide safe and convenient storage for manual and electronic pipettes using Thermo Scientific Finnpiptette Stands.



**F1 & F2 Finnpiptettes**  
Offers a variety of outstanding features including a light, smooth pipetting action. Includes add-on skin for improved grip.

[Learn More](#)

**FACTmarketing**

Repetitive-strain injuries (RSI) can occur when pipetting for extended periods of time. Weakness in your thumb or wrist when you use your pipetting hand are signs of repetitive stress injuries. Select your next pipettor wisely or ask for trial.

**Don't Forget the Tips!**  
MG Scientific offers a wide variety of universal and pipettor specific pipet tips. Tips come sterile and non-sterile, filtered and unfiltered, racked and bulk packaging.



**Finnpiptette Tips & Pipette Stand**  
1-10mL tips designed for use with Finnpiptette®, Brand®, Gilson®, Socorex®, and Labsystem® pipettors. The Finnpiptette F-stand accommodates six pipettes.

[Learn More](#)

## Hamilton - SoftGrip



**Soft Grip Fixed & Adjustable Pipets**  
New generation of pipettes with a dual indicator that improves the volume adjustment accuracy; removing the parallax errors. SoftGrip pipettes have a universal calibration key, reducing the need for costly calibration services.

[Learn More](#)



**Hamilton SoftGrip Pipette Tips**

- Made from chemically resistant, virgin polypropylene
- Safe from RNase, DNase, and pyrogens

[Learn More](#)

## Globe - Pipettors



**Single Channel & Repeating Pipets**  
The RV-Pette Pro and matching RV-Pette Pro Dispenser Syringe Tips provide flexibility and accuracy pipetting in series. New spring and seal design allows for one of the lowest plunger force requirements to help reduce fatigue.

[Learn More](#)

Search on-line to see what new and innovative products **MG Scientific** has for your lab or call **MG Customer Service** to see the many different solutions to fit your needs- 800.343.8338