



Thermo Scientific Promotions

thermo
scientific

MG Scientific, in a continued effort to keep you informed of the newest products and latest deals, would like to showcase, what our partners at Thermo Scientific™ are doing to help you get your lab up and running at peak performance. There are other promotions on new and innovative equipment to assist you and your lab as you evolve into 2022. Click here to see all current [Promotions](#), as well as updated [Covid-19 Information and Solutions](#).

Save 15% on Select Nalgene™ Consumables



1203 Griffin Low-Form Beakers

Continuously use Thermo Scientific Nalgene PMP Griffin Low-Form Beakers at temperatures of 150°C with intermittent service to 175°C. These polymethylpentene beakers have the transparency of glass, excellent chemical resistance and heat tolerance. All sizes of beakers have easy-to-read, silk-screened graduations.

1203 Griffin Low-form Beakers



5970 UnWire™ Test Tube Racks

Securely hold a full load of test tubes or centrifuge tubes for 13, 16, 20, 25 and 30mm diameter sizes with Thermo Scientific Nalgene Unwire Test Tube Racks. Molded in a single piece for exceptional strength and durability, these racks have no coating to chip, peel or rust. All sizes are molded in six bright, permanent colors. Autoclavable.

5970 UnWire Test Tube Racks



Nalgene Bottles with Closures

Excellent Chemical Resistance to most acids, bases and alcohols make these a versatile bottle. Can be used in the freezer to -100°C (-148°F) for long-term storage. Linerless caps provide the ultimate in leakproof protection. With the Polycarbonate bottles made only from high-quality, laboratory-grade plastic for dependably low leachables and extractables.

2003 Narrow-mouth LDPE Bottles

2006 Narrow-mouth PP Bottles

2103 Wide-mouth LDPE Bottles

2105 Wide-mouth PP Bottles

2110 Wide-mouth Square PP Bottles

2203 Large Narrow-mouth PP Bottles

2204 Large Amber PP Bottles

2205 Large Narrow-mouth PC Bottles

*Search on-line to see what new and innovative products **MG Scientific** has for your lab or call **MG Customer Service** to see the many different solutions to fit your needs- 800.343.8338*