

# NEW

October 2017

FALCON<sup>®</sup>


A Corning Brand

If you haven't tried **Falcon** brand pipettes yet, now is an excellent time

### Buy 3 Cases of select **Falcon** Serological Pipets, Get 1 Case Free

for new Falcon pipet customer purchases (October 15 - December 31, 2017)

- Polystyrene, disposable pipets for tissue culture, bacteriological, and research applications
- Sterile
- Noncytotoxic, nonpyrogenicity tested to less than 0.1 EU/mL
- Negative graduations for extra capacity
- Reverse graduations
- Color-coded package and markings for ease of identification
- Polyester plug to help prevent overflow
- Accurate  $\pm 2\%$  at full volume



[read more](#)

#### **P176-3 5mL Blue band**

Increments of 1/10 with an overflow capacity of 2.5mL. 200 per case

#### **P176-4 10mL Red band**

Increments of 1/10 with an overflow capacity of 3mL. 200 per case

#### **P176-5 25mL Purple band**

Increments of 0.25 with an overflow capacity of 7mL. 200 per case

#### **P176-8 25mL Purple band**

Increments of 0.5 with an overflow capacity of 11mL. 200 per case

For more Falcon-Corning products available at MG Scientific, follow this link to the [MG Scientific website](#)