



New and Improved Finnpipette F1

Adjustable and Fixed Volume
Multi-channel
GLP Kits



Improved volume adjustment

- Lighter rotation without ‘clicks’
- Improved grip on the volume adjustment button
- Updated color ring texture

Lighter plunger forces in Multichannel pipettes

- Lower spring forces on the second stop

New Appearance

- Color change to differentiate product lines
- White manifold, blue decoration ring

Improved grip

- Add-on Skin for improved grip

Improved Pipette Performance through Innovative AVG

The New and Improved F1 combines all the features of the original finnpipette with state-of-the-art innovations. The Advanced Volume Gearing (AVG) mechanism is a self-supporting modular volume adjustment mechanism that enhances pipetting performance. Separated from the pipetter body, the AVG accuracy, precision, and durability are significantly improved.



The AVG mechanism ensures a high level of accuracy and precision for the **Finnpipette F1 Multichannel Pipettes** as with the Thermo Scientific **Finnpipette F1 single-channel Adjustable-volume Pipettes** and the **single-channel Fixed-volume Pipettes**. The soft-touch tip ejection is especially useful on multichannel models, significantly reducing the effort needed to eject multiple tips. In addition, in the low volume models the super blow-out function ensures accurate dispensing for even the lowest volumes and also features an antimicrobial treatment protecting the instrument against contamination.

The **Thermo Scientific Finn timer Pipette F1 GLP Kit** comes in 4 convenient kits. Each Kit includes state-of-the-art Finn timer Pipette F1 pipettes along with everything you need for pipetting: Multiple pipettes, associated Flex tips, stand, pen, USB memory key. Click on the Specifications tab to see which Kit is right for you.

*Click on the links above to order the **NEW F1 Finn timer Pipette** that is right for your needs.*

**Inquire to Customer Service about the great deals on Discontinued F1 model Finn timer Pipettes.*