



# PROMO<sub>2019</sub>

## thermo**scientific** GENESYS— developed for users, designed for results

Built on the foundation of 60 years of spectroscopy experience, the next generation Thermo Scientific™ GENESYS™ family of Vis/UV-Vis spectrophotometers provides exceptional performance with a contemporary design.

- Start taking measurements immediately with the simple-to-use, high-resolution, color-display start menu
- Exchange samples with ease with large, accessible sample compartment
- Save your data fast any way you want it with multiple data-export options

There is a GENESYS spectrophotometer you can count on to provide accurate, reproducible results for your application requirements.



### GENESYS 30 (S284-43) Visible Spectrophotometer

Quantitative Analysis with up to Six Standards

- Better design for usability
- 6-point quant offers higher data confidence
- Scanning available with onboard software
- Backed by a 2-year warranty and skilled technical support
- Data any way you want it control

### GENESYS 40 Vis/50 (S284-71/S284-72) UV-Vis Spectrophotometers

Advanced usability, modern functionality

- 7-inch, high-resolution touchscreen interface
- Local control increases speed and reliability for routine analysis
- Large sample compartment with front access simplifies sample handling
- Wi-Fi ready— save data and methods, print results from on-board control

### GENESYS 150 (S284-74) UV-Vis Spectrophotometers

Elite performance, enhanced capabilities

- Tilttable 7-inch, high-resolution touchscreen tablet to avoid glare
- Compatible with automated cell changers, Peltier, sipper and fiber optic probe accessory options
- Wi-Fi ready— save data and methods, print results from local control

*\* Promotional pricing for GENESYS 30 (S284-43) and GENESYS 40 (S284-71) good through 6/29/2019  
GENESYS 50 (S284-72) and GENESYS 150 (S284-74) good through 3/31/2019*