



NextPette Pipette Starter Kits



The NextPette features an innovative, asymmetrical designed plunger button that can be rotated to match the comfort preference of the user. This button is freely adjustable to rotate continuously between wide setting, resulting in the more traditional feel of a round, larger diameter plunger button, or to a more narrow (thumb shaped) contour preferred by the majority of researchers.

NextPette's materials and components were carefully selected to minimize weight and ensure durability and a long service life. Internal components include precision plungers made of stainless steel (1ul, 10ul, 5000ul, and 10,000ul models) or ceramic (20ul, 200ul, and 1000ul models) and all use durable PTFE seals. The pipettors have also been designed with service in mind. They are easy to disassemble, and seals are easy to replace. Silicone lubricant is included with each pipettor. Recalibration can be done in the lab- no tools are required and the instruction manual includes full details on how to properly check pipetting accuracy and how to recalibrate if necessary.



The entire pipettor is autoclavable at 121°C and resistant to UV exposure. All pipettes are factory calibrated in accordance with ISO 8655 standards, and include a factory calibration certificate. Market leading accuracy & precision. One handed volume adjustment. Convenient volume lock. Easily visible, 4 digit display. Ergonomic, thumb shaped plunger button.



The ProPette nUVaClean™ is the first and only pipette rack that not only organizes and protects six pipettes, but also destroys up to 99% of unwanted contamination. The new nUVaClean uses a germicidal UV-C mercury vapor lamp. A high-efficiency reflector/concentrator completely bathes pipette shafts in 254nm UV light. Just place single-channel pipettes of any range or brand in the nUVaClean and press the “start” button to begin the auto-decon cycle.

Search on-line to see what standard solutions MG Scientific has for your lab or call MG Customer Service for a Special Formula that can be created to fit your needs- 800.343.8338