



SW Safety - What's Behind the Gloves?



For 35 years SW Safety has been committed to creative innovation and outstanding service; revolutionizing the healthcare and industrial products industry. Offering trusted protection inside and out, SW creates unique products that keep users safe, while making their lives easier. Innovative developments have changed the way people think, see, and feel about single use gloves.



TrueFormTM



QUALITY. COMFORT. DEXTERITY.

TrueForm is made for those who rely on daily hand protection to provide ease of use and high dexterity. With TrueForm, quality is designed into every style using exacting manufacturing processes for tip-to-cuff consistency, glove after glove.

PowerForm[®]

DESIGNED. ENGINEERED. PERFORMANCE.

PowerForm brings new hand protection technologies to workplace performance. From sweat management to enhanced grip design, PowerForm features enable workers to execute tasks with more confidence. PowerForm focuses on functional advancements to optimize the on-the-job performance every day.



MegaMan[®]



TOUGHER. ROUGHER. READY.

MegaMan pushes the boundaries of single-use hand protection performance— with new levels of abrasion resistance, sweat management, and comfort that deliver in exceedingly severe environments. When protection performance matters most, it's no wonder MegaMan gets the call for the toughest jobs.

PowerChemTM

INTRICATE. TECHNICAL. ADVANCED.

Born of SW manufacturing capabilities and advancements in material science, PowerChem offers chemical protection for highly specific applications. SW will be introducing new PowerChem styles that address a range of complex and specialized technical needs.



EcoTek

EcoTek[®] is an important technology SW has developed in support of global sustainability initiatives. SW nitrile gloves with EcoTek biodegrade in landfills through the action of naturally occurring microorganisms. According to ASTM D5526-12 biodegradation testing, biodegradation of gloves with EcoTek occur in only 945 days. Gloves break down 92.6% over 2.5 years without any impact on glove performance. MG Scientific strives to provide a more Green Initiative solution for our customers

Search on-line to see what standard solutions MG Scientific has for your lab or call MG Customer Service for a Special Formula that can be created to fit your needs- 800.343.8338