



General Purpose Refrigerators



The Thermo Scientific™ TSG Series general purpose laboratory refrigerators and chromatography refrigerators are designed with your needs in mind. The TSG series is ideal for the demanding requirements of daily use in research, clinical, and industrial laboratories storing reagents, media, proteins, and samples used daily in the lab. These TSG Series general purpose refrigerators are built with the performance, security, and quality features you have come to expect from Thermo Scientific.

TSG Series Laboratory Refrigerators

TSG Series general purpose refrigerators are designed with your laboratory in mind, combining the performance, security, and quality required for everyday laboratory applications.

- Featuring a wide range of cabinet sizes from 12-72 cu.ft.
- Whisper-quiet operation at 51 dBA
- ENERGY STAR certified Temperature, door ajar, and power failure (with battery back-up) alarms
- Green, natural refrigerants
- Adaptive defrost
- Glass (double pane) or solid door configurations
- Choice of white painted interior and exterior, or stainless steel interior and gray exterior
- LED interior lighting
- Temperature, door ajar, and power failure (with battery back-up) alarms
- Digital controller
- Low heat output
- Rear access port
- Keyed door lock(s)
- Remote alarm contacts

TSG Series Chromatography Refrigerators

TSG Series Chromatography Refrigerators are designed with your chromatography applications in mind.

- Featuring a wide range of cabinet sizes from 30-72 cu.ft.
- Hydrocarbon safety sensor ensures cargo and lab personnel are kept safe
- Whisper-quiet operation at 51 dBA
- ENERGY STAR certified
- Temperature, door ajar, and power failure (with battery back-up) alarms
- Green, natural refrigerants
- Adaptive defrost
- Double pane glass doors
- Choice of white painted interior and exterior, or stainless steel interior and gray exterior
- LED interior lighting
- Digital controller
- Low heat output
- Rear access port, plus two 2-inch side access ports
- Keyed door lock(s)
- Remote alarm contacts



Search on-line to see what standard solutions MG Scientific has for your lab or call MG Customer Service for a Special Formula that can be created to fit your needs- 800.343.8338