



# Trade-up to ENERGY STAR<sup>®</sup> certified cold storage!

More capacity, less noise and less energy than ever before!



*These TSX Series Ultra-low Freezers are “SNAP compliant”, meaning we have already transitioned from HFCs in our insulating foam to the approved natural refrigerants— well in advance of any phase-out deadlines, as we are committed to making our cold storage products support the Climate Action Plan.*



### TSX and Advanced Series Ultra-low Freezers R208-51 through R208-58

- Three sizes: 400, 600 or 700 box capacity (2-inch box)
- Energy savings: just 6.5 kWh/day (TSX400) and 7.8 kWh/day (TSX600) at -70°C setpoint
- Quiet: only 47.4 dBA (TSX400), 46.9 dBA (TSX600) and 48.0 dBA (TSX700)

### TSX Series High-performance Lab Refrigerators R209-4 through R209-11

- Quiet: just 52 dBA—up to 10x quieter than previous generation
- V-drive: significant reduction in energy with alternate compressor technology
- Ergonomic design: slanted controller to help prevent neck strain and comfortable door handle

### TSX Series High-performance Lab Freezers R208-59 through R208-66

- Two choices: -30°C auto defrost; -20°C manual defrost
- Quiet: just 52 dBA for -30°C and 50 dBA for -20°C
- V-drive: significant reduction in energy with alternate compressor technology

*All trade-up models are manufactured in a zero-waste facility and contain natural, hydrocarbon refrigerants that are compliant with the U.S. Environmental Protection Agency’s Significant New Alternatives Policy (SNAP)*