



MG Scientific - The Importance of Proper Lab Sanitation



The importance of proper lab sanitation cannot be understated. Properly cleaned labs not only protect employee safety but also create an environment conducive to producing highly accurate and repeatable results. To ensure your workspace and analytical tools are prepared to support operational efficiency, MG Scientific offers the cleaning supplies you need.

Thermo Scientific DNA Surface Decontaminant



DNA AWAY™ Refill Bottles

Suitable for gel boxes, pipets, benchtops, thermal cyclers, or other apparatus.

High-Tech Conversions Pre-saturated Wipes



FREE-SAT™ 100% Polyester

Free of additives and feature a very high level of absorbency. Individually laser cut, heat sealed-edges eliminate loose fibers.

NALGENE 5245 Pipet Washer-Rinsers



High-density Polyethylene

Wide range of cycling speeds; operates with water input from 1-½ to 12 liters per minute without "siphon-stall."



Clean-Wipes™ Universal Cleanup Wipes

Ideal for wiping optical parts, delicate glassware, electrodes, cuvettes, microscopes, or lenses. Lint-free.



Autoclavable, Reusable Safety Lab Mats

Ideal for every type of laboratory and more specialized areas. Helps facilitate workflow processes. Reversible.



Bio-Wand™ Personal UV Sanitizer

UV light (265 to 280nm) to disinfect any object/surface in the lab and beyond.

FACT miscellaneous

UV Light kills bacteria by damaging their DNA. By disrupting the chemical bonds that hold the DNA together, it prevents bacteria from being able to repair the damage, ultimately leading to the bacteria's death.



Hamilton Heated Syringe Cleaner

Do not heat the syringes with cemented needles above 50°C. Thermal expansion of the needle will fracture the glass. Do not allow the glass barrel to come into contact with the hot septum cap.



FlaskScrubber Glassware Washers

Designed to wash and dry a wide variety of narrow-neck and general purpose labware. Forced air drying through the spindles and chamber ensures thorough drying inside and out.

MG PROMOTIONS

Search on-line to see what new and innovative products MG Scientific has for your lab or call MG Customer Service to see the many different solutions to fit your needs- 800.343.8338